**EN-1**: Single use nitrile examination glove. Powder Free. Not Made With Natural Rubber Latex. Textured Finish. Non-Sterile. Contains Potential Type I Allergens. Protection is limited to the hand only. If irritation occurs seek medical advice. **Disposal:** Once discarded no re-use. **Storage:** Once opened no re-use. **Use:** Within 3 years from date of manufacture. Service life cannot be specified and depends on application and responsibility of user to ascertain suitability of the glove for its intended use. Before usage, inspect the gloves for any defect or impregnations. **Inappropriate glove use:** Tactile sensitivity and dexterity reduces sharply with a poor fitting glove, which can cause fatigue in the hands and fingers. Wrong size and poor fitting lead to poor hand protection. **Glove donning:** Choose the appropriate glove size. Check for any physical damage or contamination prior to each use of the glove. Discard if there is any damage or contaminant. **Glove doffing:** Remove the glove as soon as you see wear or damage. When the glove is covered in a contaminant that is not removable or presents a potential hazard, it should be left as is and right gloves off alternately using the gloved hand so that the gloves are removed without the contaminant contacting the skin. **Hand Hygiene:** To maintain hand hygiene wash hands thoroughly before donning gloves and after removing gloves. **Gloves:** Gloves should be disposed by incineration or via the regulation of authorities. To access and download Declaration of Conformity please visit [g-drc.com/skylite](http://g-drc.com/skylite). The product is compliant with the PPE Regulation EU 2016/425, certified according to EN 12941, A2012/A2013.

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Population levels indicate the change in Percentage Deviations of the class after exposure to the challenge chemical.

27

16	Performance Level	AQL	Inspection Levels	Tested in accordance to EN ISO 374-2:2014 Penetration air leak: PASS Penetration water leak: PASS Performance level: Level 2
	Level 3	< 0.65	GI	Micro-organisms risk in accordance with EN ISO 374-5:2016
	Level 2	< 1.5	GI	Protection against bacteria and fungi: PASS
	Level 1	< 4.0	S4	Protection against viruses: PASS

Tahrif Intended: Medical and RPE Cat II

Comply to MDR2017/745 Class I

**STATEMENT AND CAUTION** Test results are taken from palm area of the glove. Glove performance quoted is based on laboratory data and may not reflect the actual duration of protection in the workplace due to the differentiation between mixtures and pure chemicals and other factors influencing the performance such as temperature, abrasion, degradation etc. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. When used, protective gloves may not give less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. There are no possible allergies within the glove that are known to cause harm to the wearer.

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**LOT**

**EU REP** Globus EMEA Ltd Authorized Representative in the European Commun

Journal of Oral Rehabilitation 2006 33: 103–109

CLORUS

Manufactured for:  
**Globus (Shetland) Ltd.**

12 Trafford Point, Twining Road, Trafford Park, Manchester M17 1SH, UK | T: +44 (0)161 877 4747 | [www.gldbusgroup.co](http://www.gldbusgroup.co)

Manufacturer: GLOBUS EMEA FZE, Jafza View 1, Office 220, Jebel Ali, PO Box 61195, Dubai, United Arab Emirates