

Haika HKC02 Surgical Gown

Reusable PVC Coated Surgical Gown

Key Features

Ergonomic fit for freedom of movement

Secure back neck studded adjustment

Elasticated soft touch ribbed cuffs

Thumb loops for added security

One size fits all wraparound design for enhanced coverage

Easily adjusted with side tie fastening

Fluid repellent PVC coating with soft comfortable feel

Machine washable up to 75°C for 100 washes 'under testing'

Lighter, more flexible than standard PVC

Stitched and high frequency top welded seams for greater protection

Quick drying



Made in UK



Medical device

EN13795-1

Ordering Information

- **Product Code:** HKC00002EZ
- **Pack Quantity:** 1 singles per bag
- **Case Quantity:** 10 singles per carton

EU Compliance

- Class 1 Medical Device
- EN13795-1 Standard Performance

Haika HKC02 Surgical Gown

STORAGE AND MAINTENANCE

WASH INSTRUCTIONS

- Non-ionic detergent should be used
- Do not use bleach
- Do not tumble dry
- Do not iron
- Do not dry clean
- Wash at a maximum temperature of 75°C
- Hang dry to minimise creasing
- Maximum 100 washes (Under testing)

STORAGE

- It is recommended that the garments are stored in their original packaging clear of UV sources. Following use and cleaning, garments should be stored in a cool, dry and clean environment

Components

- **Body:** 81% PVC 19% Polyester
- **Cuff:** Polyester
- **Thumb Loop:** Braided Elastic
- **Tie:** Polyester Herringbone Weave
- **Binding:** Polyester Blend

Sizing

- **Item:** One size
- **Center Back Length (cm):** 120
- **Width (cm):** 160
- **Under Arm Length including cuff (cm):** 64.5

More information:

EUROPE

www.globusgroup.com

E: sales@globus.co.uk
T: +44 (0)161 877 4747
F: +44 (0)161 877 4746

MIDDLE EAST AND AFRICA

www.globusgroup.com/gcc

E: gcc@globusgroup.com
T: +971 4 882 9962
F: +971 4 882 9963

